

“DIAQUICK” Multi-Drug Panels

for human urine samples

Multi-3 Drug Panel	BZO,COC,MOP
- REF Z06576CE	Cont.: 30 panels, individually packed (30x REF Z06576B)
- REF Z06576B	Cont.: 1 panel, individually packed
Multi-3/1 Drug Panel	BUP, MOP, MTD
- REF Z09577CE	Cont.: 30 panels, individually packed (30x REF Z09577B)
- REF Z09577B	Cont.: 1 panel, individually packed
Multi-4 Drug Panel	AMP,COC,MOP,THC
- REF Z02575CE	Cont.: 30 panels, individually packed (30x REF Z02575B)
- REF Z02575B	Cont.: 1 panel, individually packed
Multi-5 Drug Panel	BZO,COC,MET,MOP,THC
- REF Z05236CE	Cont.: 30 panels, individually packed (30x REF Z05236B)
- REF Z05236B	Cont.: 1 panel, individually packed
Multi-5/3 Drug Panel	AMP,COC,MET,MOP,THC
- REF Z06502CE	Cont.: 30 panels, individually packed (30x REF Z06502B)
- REF Z06502B	Cont.: 1 panel, individually packed
Multi-5/4 Drug Panel	AMP,COC,MDMA,MOP,THC
- REF Z11504CE	Cont.: 30 panels, individually packed (30x REF Z11504B)
- REF Z11504B	Cont.: 1 panel, individually packed
Multi-5/6 Drug Panel	AMP,BZO,COC,MOP,THC
- REF Z06506CE	Cont.: 30 panels, individually packed (30x REF Z06506B)
- REF Z06506B	Cont.: 1 panel, individually packed
Multi-6 Drug Panel	BZO,COC,MET,MOP,MTD,THC
- REF Z98907CE	Cont.: 30 panels, individually packed (30x REF Z98907B)
- REF Z98907B	Cont.: 1 panel, individually packed
Multi-6/1 Drug Panel	AMP,BZO,COC,MET,MOP,THC
- REF Z03220CE	Cont.: 30 panels, individually packed (30x REF Z03220B)
- REF Z03220B	Cont.: 1 panel, individually packed
Multi-6/3 Drug Panel	BUP, BZO, COC, MTD, OPI, THC
- REF Z08930CE	Cont.: 30 panels, individually packed (30x REF Z08930B)
- REF Z08930B	Cont.: 1 panel, individually packed
Multi-6/4 Drug Panel	AMP,BUP,BZO,MET,MOP,THC
- REF Z08940CE	Cont.: 30 panels, individually packed (30x REF Z08940B)
- REF Z08940B	Cont.: 1 panel, individually packed
Multi-6/6 Drug Panel	BUP,COC,MET,MOP,MTD,THC
- REF Z13960CE	Cont.: 30 panels, individually packed (30x REF Z13960B)
- REF Z13960B	Cont.: 1 panel, individually packed
Multi-6/7 Drug Panel	BUP,BZO,COC,MOP,MTD,THC
- REF Z09970CE	Cont.: 30 panels, individually packed (30x REF Z09970B)
- REF Z09970B	Cont.: 1 panel, individually packed
Multi-6/10 Drug Panel	AMP,BZO,COC,MOP,MTD,THC
- REF Z11911CE	Cont.: 30 panels, individually packed (30x REF Z11911B)
- REF Z11911B	Cont.: 1 panel, individually packed
Multi-7 Drug Panel	AMP,BUP,BZO,COC,MTD,MOP,THC
- REF Z12730CE	Cont.: 30 panels, individually packed (30x REF Z12730B)
- REF Z12730B	Cont.: 1 panel, individually packed
Multi-10 Drug Panel	AMP,BAR,BZO,COC,MDMA,MET,MOP,MTD,TCA,THC
- REF Z04230CE	Cont.: 30 panels, individually packed (30x REF Z04230B)
- REF Z04231CE	Cont.: 10 panels, individually packed (10x REF Z04230B)
- REF Z04230B	Cont.: 1 panel, individually packed
Multi-10/1 Drug Panel	AMP,BAR,BZO,BUP,COC,MDMA,MET,MOP,MTD,THC
- REF Z06235CE	Cont.: 30 panels, individually packed (30x REF Z06235B)
- REF Z06236CE	Cont.: 10 panels, individually packed (10x REF Z06235B)
- REF Z06235B	Cont.: 1 panel, individually packed
Multi-10/2 Drug Panel	AMP,BAR,BZO,COC,MDMA,MOP,MTD,OPI,TCA,THC
- REF Z06102CE	Cont.: 30 panels, individually packed (30x REF Z06102B)
- REF Z06102B	Cont.: 1 panel, individually packed
Multi-10/3 Drug Panel	AMP, BZO, COC, MDMA, MOP, MTD, OPI, PCP, TCA, THC
- REF Z06103CE	Cont.: 30 panels, individually packed (30x REF Z06103B)
- REF Z06103B	Cont.: 1 panel, individually packed
Multi-10/4 Drug Panel	AMP,BAR,BUP,BZO,COC,MDMA,MET,MTD,OPI,THC
- REF Z06104CE	Cont.: 30 panels, individually packed (30x REF Z06104B)
- REF Z06104B	Cont.: 1 panel, individually packed
Multi-10/5 Drug Panel	AMP,BAR,BZO,BUP,COC,MET,MOP,MTD,TCA,THC
- REF Z06105CE	Cont.: 30 panels, individually packed (30x REF Z06105B)
- REF Z06105B	Cont.: 1 panel, individually packed
Multi-10/6 Drug Panel	AMP,BAR,BZO,COC,MET,MOP,MTD,PCP,TCA,THC
- REF Z06106CE	Cont.: 30 panels, individually packed (30x REF Z06106B)
- REF Z06106B	Cont.: 1 panel, individually packed
Multi-10/7 Drug Panel	AMP,BAR,BZO,COC,MET,MOP,MTD,PCP,TCA,THC
- REF Z06107CE	Cont.: 30 panels, individually packed (30x REF Z06107B)
- REF Z06107B	Cont.: 1 panel, individually packed

All products contain a package insert!

For in vitro diagnostic use only. For use by medical professionals only.
 For diagnosis and therapeutic monitoring only.

INTENDED USE

The “DIAQUICK” Multi-Drug Panels (urine) are rapid, lateral flow chromatographic immunoassays for the simultaneous, qualitative detection of the following drugs and their metabolites:

Parameter	Short	Calibrator Substance	Cut-off
Amphetamine	AMP	d-Amphetamine	1,000 ng/mL
Barbiturates	BAR	Secobarbital	300 ng/mL
Buprenorphine	BUP	Buprenorphine	10 ng/mL
Benzodiazepines	BZO	Oxazepam	300 ng/mL
Cocaine	COC	Benzoyllecgonine	300 ng/mL
Ecstasy	MDMA	d,l-Methylenedioxymethamphetamine	500 ng/mL
Methamphetamine	MET	d-Methamphetamine	1,000 ng/mL
Methodone	MTD	Methodone	300 ng/mL
Opiate, Morphine, Heroine	MOP	Morphine	300 ng/mL
Opiate, Morphine, Heroine	OPI	Morphine	2,000 ng/mL
Phencyclidine	PCP	Phencyclidine	25 ng/mL
Tricyclic Antidepressants	TCA	Nortriptyline	1,000 ng/mL
Marihuana/Cannabis	THC	11-nor- Δ^9 -THC-9-COOH	50 ng/mL

This test will detect other related compounds, please refer to the Analytical Specificity table in this insert. This assay provides only a preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are obtained. For in vitro diagnostic use only

TEST PRINCIPLE

The “DIAQUICK” Multi-Drug Panels (urine) are immunoassays based on the principle of competitive binding. Drugs which may be present in the urine specimen compete against their respective drug conjugate for binding sites on their specific antibody. During testing, a urine specimen migrates upward by capillary action. A drug, if present in the urine specimen below its cut-off concentration, will not saturate the binding sites of its specific antibody coated on the particles. The antibody coated particles will then be captured by the immobilized drug conjugate and a visible colored line will show up in the test line region of the specific drug strip. The colored line will not form in the test line region if the drug level is above its cut-off concentration because it will saturate all the binding sites of the antibody coated on the particles. A drug-positive urine specimen will not generate a colored line in the specific test line region of the strip because of drug competition, while a drug-negative urine specimen or a specimen containing a drug concentration less than the cut-off will generate a line in the test line region. To serve as a procedural control, a colored line will always appear at the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred.

WARNINGS AND PRECAUTIONS

- For medical and other in vitro diagnostic use only. Do not use after the expiration date.
- The test panel should remain in the sealed pouch until use.
- All specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The used test panels should be discarded according to federal, state and local regulations.

REAGENTS

Each test line contains anti-drug mouse monoclonal antibody and corresponding drug-protein conjugates. The control line contains goat anti-rabbit IgG polyclonal antibodies and rabbit IgG.

STORAGE

The “DIAQUICK” Multi-Drug Panels can be stored refrigerated or at room temperature (2-30°C). The test panel is stable through the expiration date printed on the sealed pouch. The test panel must remain in the sealed pouch until use. DO NOT FREEZE. Do not use beyond the expiration date.

SAMPLE COLLECTION AND PREPARATION


The urine must be collected in a clean and dry container. Urine collected at any time of the day may be used. Urine specimens exhibiting visible precipitations should be centrifuged, filtered or allowed to settle to obtain a clear specimen for testing. Urine specimens may be stored at 2-8°C for up to 48 h prior to testing. For prolonged storage, specimens may be frozen and stored below -20°C. Frozen specimens should be thawed and mixed before testing.

ASSAY PROCEDURE

Allow the test panel, urine specimen, and/or controls to equilibrate to room temperature (15-30°C) prior to testing

1. Remove the test panel from the sealed pouch and use it as soon as possible.

2. Take off the protective cap plugged on the test panel. With arrows pointing towards the urine specimen, immerse the test panel vertically into the urine specimen for 10-15 seconds. Do not allow the urine sample to touch the plastic cassette when immersing the test device into the urine sample. Avoid immersion of the cassette deeper than the mark indicated with the arrows on the device and avoid any direct contact of the sample with the test region.



3. Put the protective cap back onto the test panel. Place the test panel on a non-absorbent flat surface, start the timer and wait for the red line(s) to appear. Read the results at 5 minutes. Do not interpret results after 10 minutes.

INTERPRETATION OF RESULTS

NEGATIVE: A colored line in the control region (C) and a colored line in the test line region (T) for a specific drug indicate a negative results. This indicates that the drug concentration in the urine specimen is below the designated cut-off level for that specific drug.

*NOTE: The shade of color in the test region (T) may vary, but it should be considered negative whenever there is even a faint pink line.

POSITIVE: A colored line in the control line region (C) but no line in the test line region (T) for a specific drug indicates a positive results. This indicates that the drug concentration in the urine specimen exceeds the designated cut-off level.

INVALID: Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test using a new test panel. If the problem persists, discontinue using the lot immediately and contact your local distributor

QUALITY CONTROL

A procedural control is included in the test. A red line appearing in the control region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique. Control standards are not supplied with this kit; however, it is recommended that positive and negative controls be tested as good laboratory practice to confirm the test procedure and to verify proper test performance.

LIMITATIONS

1. The “DIAQUICK” Multi-Drug Panels (urine) provide only a preliminary analytical result. A more specific chemical method must be used to obtain a confirmed result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method.
2. It is possible that technical or procedural errors, as well as other interfering substances in the urine specimen may cause erroneous results.
3. Adulterants, such as bleaching agents in urine specimens may produce erroneous results regardless of the analytical method used. If adulteration is suspected, the test should be repeated with another urine specimen.
4. A positive result indicates presence of the drug or its metabolites but does not indicate

